Instruction  
**WIESLAB® Complement system**  
Classical pathway  
Qualitative and Semi-Quantitative test

For in vitro diagnostic use

Enzyme immunoassay for assessment of  
Complement functional activity

- Break apart microtitration strips (12x8) 96 wells
- Store the kit at +2-8° C
- Store the positive control and activity control at -20° C

Doc. no. LABEL-DOC-0031 7.0

Effective Date: 22-Jan-2026

**WIESLAB® Complement system**  
**Classical pathway**

English: page ..... 1

[Czech: page ..... 17](#)

[Estonian: page ..... 33](#)

[French: page ..... 50](#)

[German: page ..... 67](#)

[Hungarian: page ..... 85](#)

[Icelandic: page ..... 103](#)

[Italian: page ..... 119](#)

[Norwegian: page ..... 137](#)

[Portuguese: page ..... 153](#)

[Serbian: page ..... 171](#)

[Spanish: page ..... 188](#)

[Swedish: page ..... 206](#)

## 1. INTENDED USE

The WIESLAB® Complement system Classical pathway (CP) kit is an enzyme immunoassay for the qualitative and/or semi-quantitative determination of functional classical complement pathway in human serum. The analysis should be performed by trained laboratory professionals.


FOR IN VITRO DIAGNOSTIC USE.

## 2. SUMMARY AND EXPLANATION

The complement system plays an essential role in chronic, autoimmune and infectious disease. There are three pathways of complement activation (Figure 1), namely the classical, the MBL (lectin) and the alternative pathway.

Impaired complement activity causes humans to become susceptible to repetitive fulminant or severe infections and may contribute to development of autoimmune disease. Inappropriate activation of complement contributes to chronic inflammation and tissue injury.

*In vitro* activation of the complement cascade leads to the consumption of complement components which, in turn, leads to a decrease in their concentration. Thus, the determination of complement proteins or complement activity is used to indicate whether the complement system has been activated by an immunologic and/or pathogenic mechanism. Both functional and immunochemical complement measurements are used to evaluate patients when a complement-activating disease is suspected or an inherited deficiency is possible. The level of complement activity evaluated by functional assays, such as *WIESLAB* Complement kits, take into account the rate of synthesis, degradation, and consumption of the components and provides a measure of the integrity of the pathways as opposed to immunochemical methods which specifically measure the concentration of various complement components.

**Figure 1:** Schematic figure of the complement system. The complement system can be activated via the classical, MBL (Lectin), or alternative pathway by different activating molecules. Upon activation C3 convertases are formed which in turn cleave C3 and subsequently C5 can be cleaved. Cleavage of C5 to C5b initiates the formation of the MAC complex (also called terminal complement complex (TCC) or C5b-9).

### 3. PRINCIPLE OF WIESLAB® COMPLEMENT SYSTEM CLASSICAL PATHWAY

The WIESLAB® Complement system Classical pathway assay combines principles of the hemolytic assay for complement activation with the use of labelled antibodies specific for a neoantigen produced as a result of complement activation. The amount of neoantigen generated is proportional to the functional activity of the classical pathway.

In the Complement system Classical pathway kit, the wells of the microtitre strips are coated with a specific activator of the classical pathway. The plate in combination with sample dilution buffer composition and patient serum dilution level ensure that only the classical pathway is activated. During the incubation of the diluted patient serum in the wells, complement is activated by the specific coating.

The wells are then washed and the amount of C5b-9 complex formed on the plate surface is detected with a specific alkaline phosphatase labelled antibody to the C5b-9 neoantigen formed during Membrane Attack Complex (MAC) formation.

After a further washing step, detection of specific antibodies is obtained by incubation with alkaline phosphatase substrate solution. The amount of complement activation correlates with the colour intensity and is measured in terms of absorbance (optical density, OD).

#### Metrological traceability

The value assigned to the positive control is defined as 100 % of a normal complement pathway activity and is not traceable to SI units. It is assigned using pooled human normal sera assessed with the described ELISA assay procedure with an internal reference system as the highest calibration hierarchy level.

### 4. WARNINGS, CAUTIONS AND PRECAUTIONS

- FOR *IN VITRO* DIAGNOSTIC USE.
- The human serum components used in the preparation of the controls in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come in to contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- Safety data sheet for all hazardous components contained in this kit is available on request from Svar Life Science.
- Do not use components past the expiration date.
- Kit reagents are kit specific and should not be mixed between kit lots.
- Dispose of all used components, used/leftover samples and leftover controls as biohazardous waste according to local regulations. Dispose of other leftover components as hazardous waste according to local regulations.
- Any serious incident that has occurred in relation to this device shall be reported to Svar Life Science AB and the competent authority of the EU Member State or country in which the user/patient is established.

Wash Solution (30x Conc.)

#### WARNING

**Contains:** Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)

H317 May cause an allergic skin reaction.  
H412 Harmful to aquatic life with long lasting effects.  
P261 Avoid breathing spray.  
P273 Avoid release to the environment.  
P280 Wear protective gloves.  
P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

Dilution Buffer CP, Conjugate Solution, Negative Control, Substrate pNPP

EUH208 Contains "Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)" May produce an allergic reaction.

EUH210: Safety data sheet available on request.

## 5. KIT CONTENTS

- 1 sealed microtitre plate with break-apart wells (12x8) coated with human IgM.
- 2 vials (35 mL) Diluent CP (Dil CP), labelled blue.
- 1 vial (0.2 mL) negative control (NC) containing human serum (to be diluted as for a patient serum sample).
- 1 vial (0.2 mL) positive control (PC) containing lyophilized human serum, see section 9.3 "Reconstitution of positive control", below.
- 1 vial (vol, see CoA) activity control (AC) for semi-quantitative application, containing lyophilized human serum (different origin than PC), see section 9.3.1 "Reconstitution of activity control" under procedure for semi-quantitative application.
- 1 vial (13 mL) conjugate containing alkaline phosphatase-labelled antibodies to C5b-9 (blue colour).
- 1 vial (13 mL) Substrate solution pNPP, ready to use.
- 1 vial (30 mL) wash solution, 30x concentrated.

**Please note:** The reconstitution volume for AC is indicated in Certificate of Analysis (CoA) (XXX µl) and on AC label.

## 6. MATERIALS OR EQUIPMENT REQUIRED BUT NOT PROVIDED

- Microplate reader with filter 405 nm.
- Precision pipettes with disposable tips.
- Washer for strips, absorbent tissue, tubes and a timer.
- Incubator capable of maintaining 37°C. If a CO<sub>2</sub> cabinet is used, make sure that the CO<sub>2</sub> supply is disconnected/off

## 7. HANDLING AND STORAGE

- The reagents should be stored at 2-8°C except the positive and the activity control.
- All reagents in the kit are ready to use except wash solution and controls.
- The positive and activity controls should be stored at -20°C until reconstitution.
- The reconstituted positive and activity control should be stored at < -70°C and may be thawed once.

## 8. SAMPLE PREPARATION

This test is performed on serum specimens. Blood samples are to be collected using aseptic venipuncture technique and serum obtained using standard procedures. A minimum of 5 mL of whole blood is recommended. Allow blood to clot in serum tubes, for 60-65 minutes at room temperature (20-25°C). Centrifuge blood samples and transfer cell-free serum to a clean tube.

Sera must be properly handled to prevent *in vitro* complement activation. Sera should be frozen at -70°C or lower in tightly sealed tubes for extended storage or for transport on dry ice. Samples should not be frozen and thawed more than once.

Do not use sera which are icteric, lipemic and haemolysed. Heat-inactivated sera cannot be used. Plasma cannot be used. The Clinical and Laboratory Standards Institute provides recommendations for storing blood specimens (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

## 9. PREPARATION AND HANDLING OF REAGENTS

### 9.1 Preparation of washing solution

In case salt crystals are observed in the vial with concentrated wash solution, place the vial at 37°C in a water bath until the crystals have dissolved before dilution of wash solution.

Dilute 30 mL of the 30x concentrated wash solution in 870 mL distilled water. When stored at 2-8°C, the diluted wash solution is stable until the date of expiration of the kit.

### 9.2 Negative Control (NC)

NC should be diluted as a patient serum sample.

### 9.3 Reconstitution of positive control (PC)

Reconstitute the PC according to the following procedure:


1. Gently tap down all lyophilized material to the bottom of the vial and remove the cap.
2. Immediately add 200 µL of distilled water directly to the lyophilized material.
3. Replace the cap.
4. Allow the vial to stand on ice for 5 minutes and then gently shake or vortex occasionally until completely dissolved.
5. For further dilution, See section 9.3.1 or 9.3.2 for semi-quantitative and qualitative protocol, respectively.

The reconstituted PC can be stored for up to 4 hours prior to use if kept at 2-8°C or on ice. It should be stored at <-70°C and may be thawed once.

#### 9.3.1. Semi-quantitative protocol

##### Preparation of calibration curve

For dilution of reconstituted PC to calibrators please see Figure 2 and Table 1. The calibrator can be left at room temperature up to 1 h before use. The calibrator must be prepared fresh and cannot be stored for later usage.

**Figure 2:** Dilution scheme for calibrator curve. Reconstituted positive control (PC) is mixed with kit diluent (DIL) according to the table.

### Reconstitution of activity control (AC)

1. Gently tap down all lyophilized material to the bottom of the vial and remove the cap.
2. Add immediately the volume distilled water indicated in the CoA /AC label directly to the lyophilized material.
3. Reattach the cap.
4. Allow the vial to stay on ice for 5 minutes and then shake or vortex gently until complete dissolution.
5. Dilute the reconstituted control in the same way as a patient serum sample.

The reconstituted activity control can be stored for up to 4 hours prior to use if kept at 2-8°C or on ice. It should be stored at <-70°C and may be thawed once.

### 9.3.2. Qualitative protocol

Dilute the reconstituted positive control according to the instructions in section 9.4: Dilution of samples.

### 9.4 Handling of Serum

Partially thaw frozen sera by briefly placing it in a 37°C water bath with gentle mixing. After partially thawing, immediately place the tubes on ice and leave there until completely thawed. Mix briefly on a vortex mixer.

### Dilution of samples

Dilute the serum 1/101 with Diluent CP, blue label, (500 µL Diluent + 5 µL serum) and mix thoroughly but gently on a vortex. The diluted serum can be left at room temperature for a maximum of 60 minutes before analysis.

## 10. ASSAY PROCEDURE

### 10.1. Rinsing protocol

Empty the wells and wash 3 times with 300  $\mu$ L washing solution, filling and emptying the wells each time. After the last wash, empty the wells by tapping the strip on an absorbent tissue.

### 10.2 Plate layout

#### 10.2.1 Semi-quantitative protocol

Pipet 100  $\mu$ L/well in duplicate of calibrator (100%-0%), NC, AC and diluted serum samples (P) according to Table 2. The 0% calibrator should be used as a blank.

**Table 2:** Suggested plate layout, semi-quantitative protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	100%	12.5%	P1									
B	100%	12.5%	P1									
C	75%	0%	P2									
D	75%	0%	P2									
E	50%	NC	etc									
F	50%	NC										
G	25%	AC										
H	25%	AC										

#### 10.2.2. Qualitative protocol

Pipet 100  $\mu$ L/well in duplicate of Diluent as a blank, PC, NC and diluted serum samples (P) according to Table 3.

**Table 3:** Suggested plate layout, qualitative protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	Dil CP	P2										
B	Dil CP	P2										
C	PC	etc										
D	PC											
E	NC											
F	NC											
G	P1											
H	P1											

### 10.3. Assay protocol

Remove only the number of wells needed for testing, resealing the aluminium package carefully. Let all solutions equilibrate to room temperature (20-25°C) before analysis.

#### 1. Incubation of samples

Fill the plate according to section 10.2. Incubate for 60-70 minutes at +37°C with lid. Please note that no incubation should be performed under CO<sub>2</sub> atmosphere. If a CO<sub>2</sub> cabinet is used, make sure that the CO<sub>2</sub> supply is disconnected/off.

#### 2. Wash

Wash the plate according to section 10.1: Rinsing protocol.

#### 3. Conjugate incubation

Add 100 µL conjugate to each well. Incubate for 30 minutes at room temperature (+20-25°C).

#### 4. Wash

Wash 3 times as before.

#### 5. Substrate incubation

Add 100 µL substrate solution to each well, incubate for 30 minutes at room temperature (+20-25°C).

#### 6. Reading the plate

Read the absorbance (OD) at 405 nm on a microplate reader.

Optional: 5 mM EDTA can be used as stop solution, 100 µL/well. Read the absorbance of the wells within 60 minutes.

## 11. QUALITY CONTROL


Certificate of Analysis (CoA) included in the kit is lot specific and is to be used to verify results obtained by our laboratory. The results indicated on CoA are to be used as a guideline only. The results obtained by your laboratory may differ.

The OD of the PC should be >1.0 and the OD of the NC <0.2 after subtraction of the blank (Dil CP, 0%). If any of the controls are not within their respective ranges, the test should be considered invalid and then repeated.

## 12. INTERPRETATION OF RESULTS

### 12.1. Semi-quantitative protocol

Subtract the absorbance of the 0 % - calibrator from all OD values. Curve fit 4-parameter logistic (Marquardt) is recommended. See an example of a standard curve in Figure 3.

**Figure 3:** Example of a standard curve. The figure above shows an example of a semi-quantitative standard curve and should not be used for actual patient sample interpretation.

In cases where the obtained sample values are higher than the highest Calibrator 100%, the samples can be diluted 1/201 and retested. Please note that the obtained activity value in this case should be adjusted according to applied sample dilution.

It is recommended that each laboratory establishes its own reference level and cut-off value for deficiencies.

## 12.2. Qualitative protocol

It is recommended to use the OD results to calculate complement activity (% of PC) values for ease of interpretation. This calculation also facilitates comparisons between different analyses as OD is more prone to change from run to run. The method used to calculate complement activity is described below.

1. Subtract the absorbance of the blank (Diluent) from the absorbances of the NC, PC and the samples.
2. Calculate the mean OD values for the NC, PC and samples.
3. Calculate complement activity using the formula below:

$$\text{Complement activity (\% of PC)} = \frac{\text{Sample} - \text{NC}}{\text{PC} - \text{NC}} \times 100$$

The NC and PC are intended to monitor for substantial reagent failure. The PC will not ensure precision at the assay cut-off. It is recommended that each laboratory establishes its own reference level and cut-off value for deficiencies.

A negative result, i.e. deficiency, should always be verified by testing a new sample to ensure that no inadvertent *in vitro* complement activation has taken place prior to execution of the assay.

## 13. EVALUATION CRITERIA

### 13.1. Expected Results

When decreased levels of complement components or complement function are found, a deficiency or an ongoing, immunologic process, leading to increased breakdown of components and depression of complement levels is considered by clinicians. Increased complement levels are usually a nonspecific expression of an acute phase response.

The WIESLAB® Complement system Classical pathway can be helpful for detection of complement deficiencies related to the classical pathway. A more complete and in-depth functional assessment of all three complement pathways may be achieved using WIESLAB® Complement system Screen kit as shown in Table 4.

**Table 4:** Combination of results in the different pathways and possible deficiencies

Classical pathway	MBL pathway	Alternative pathway	Possible deficiency
Positive	Positive	Positive	None
Negative	Positive	Positive	C1q, C1r, C1s
Positive	Positive	Negative	Factor B, D and P
Positive	Negative	Positive	MBL, MASP-2
Negative	Negative	Negative	C3, C5, C6, C7, C8, C9
Negative	Negative	Positive	C4, C2 or combination

## 13.2. Limitations

The individual patient's complement level cannot be used as a measure of disease severity, as it may vary from patient to patient. Thus, it is difficult to obtain an absolute standardisation of results.

The test should not be relied upon as the sole basis of decisions on clinical therapy but should be used in combination with clinical symptoms and the results of other available tests. Therapy should not be started on basis of the complement assay result. Initiation or changes in treatment should not be based on changes in complement levels alone, but rather on careful clinical observation.

## 14. PERFORMANCE CHARACTERISTICS

### 14.1. Clinical performance


The normal distribution within 2SD (Standard Deviation) is indicated in Figure 4 and Table 5. Results within this range indicate a normal functionality of the classical pathway. It is recommended that each laboratory confirms or establishes their own reference range for the population they serve.

A value below the reference ranges indicates either increased activation, resulting in consumption of the classical complement pathway capacity, or a genetically determined low activity.

Values below 5% strongly suggest a complete deficiency either caused by excessive activation or an inherited deficiency in the respective pathway. To establish which complement factor(s) causes the lowered activity further analysis of complement proteins is needed.

A negative result, i.e. suspected deficiency, should always be verified by testing a new, carefully handled sample to ensure that no *in vitro* complement activation has taken place.

Sera from 120 apparently healthy blood donors were tested in the WIESLAB® Complement system Classical pathway, qualitative application. For the semi-quantitative application, the normal range is based on mathematical theoretical calculation. The complement activity of the 120 apparently healthy blood donors is summarized in Figure 4 and Table 5. In the study no blood donor was below 40 %.

**Figure 4:** Histogram of reference population analysed in classical pathway, semi-quantitative (A) and qualitative (B) protocol, respectively.

**Table 5:** Values for the 120 blood donors analysed with semi-quantitative and qualitative protocol, respectively.

Protocol	n	Mean (%)	±2SD (%)	Median (%)
<b>Semi-quantitative</b>	120	89	66-113*	92
<b>Qualitative</b>	120	99	69-129*	100

\*) This is a statistical calculation and will not guarantee a true cut-off. It's recommended that each laboratory establishes its own reference level and cut-off value for suspected deficiency.

Sera with known complement deficiencies and specific complement factor depleted sera were tested in the assay (Table 6 and 7). All deficient/depleted sera were low in the assay and gave values below 15% and 5%\*\* in the semi-quantitative and qualitative protocol, respectively.

\*\* See "M.A. Seelen et. al, Functional analysis of the classical, alternative and MBL pathways of the complement system: standardization and validation of a simple ELISA. J Immunol Meth 2005, 296, 187-198", for extended tests of deficient patient samples tested with qualitative application

**Table 6:** Samples with known deficiencies

Deficiency	C2	C3	C4	C5	C7	C8
<b>Number of patients</b>	5	1	1	1	2	2
<b>Number of deficient sera detected</b>	5	1	1	1	2	2

**Table 7:** Samples with depleted components

Depletion	C1q	C3	C4	C5	C7
<b>Number of depleted sera</b>	2	1	1	1	1
<b>Number of detected depletion</b>	2	1	1	1	1

## 14.2. Precision

### 14.2.1 Inter-assay precision

#### Semi-quantitative protocol

**Table 8:** Inter-assay precision for the semi-quantitative application was determined by testing seven samples in eight replicates on three different occasions.

	1	2	3	4	5	6	7
<b>Mean value %</b>	71	73	69	72	25	35	31
<b>SD</b>	9	9	6	11	1	2	2
<b>CV %</b>	13%	13%	9%	15%	4%	5%	5%

#### Qualitative protocol

**Table 9:** Inter-assay precision for the qualitative application was determined by testing three samples in duplicate. Results were obtained for six different runs.

	1	2	3
<b>Mean value %</b>	98	92	21
<b>SD</b>	4.3	3.9	1.7
<b>CV %</b>	4	4	8

### 14.2.2 Intra-assay precision

#### Semi-quantitative protocol

**Table 10:** Intra-assay precision for the semi-quantitative application was determined by testing seven different samples in eight replicates on one occasion.

	1	2	3	4	5	6	7
<b>Mean value %</b>	79	76	72	83	24	34	30
<b>SD</b>	10	11	6	9	0	1	1
<b>CV %</b>	13%	15%	8%	10%	2%	4%	2%

#### Qualitative protocol

**Table 11:** Intra-assay precision for the qualitative application was determined by testing one sample in 40 wells.

Assay	Mean value %	SD	CV %
<b>CP</b>	85	2.9	3

### 14.2.3 Batch-to-batch variation

#### Semi-quantitative protocol

**Table 12:** Batch to batch variation in the semi-quantitative application was determined by testing seven samples in duplicate on three different batches by three different persons.

	1	2	3	4	5	6	7
<b>Mean value %</b>	8	73	85	24	37	74	78
<b>SD</b>	0.80	14.55	12.50	1.32	1.91	7.46	6.08
<b>CV %</b>	10%	20%	15%	5%	5%	10%	8%

### 14.3. Linearity

#### Semi-quantitative protocol

**Table 13:** Dilution recovery was determined by testing five serial dilutions for three different samples.

Sample	Dilution	Mean Measured Activity (%)	Theoretical Activity (%)	Dilution Corrected % Recovery
1	1/200	55	55	100
	1/400	32	28	114
	1/800	17	14	121
	1/1600	0	7	0
2	1/200	46	46	100
	1/400	25	23	109
	1/800	13	12	108
	1/1600	6	6	100
3	1/100	84	84	100
	1/200	37	42	88
	1/400	21	21	100
	1/800	11	11	100
	1/1600	7	6	117

### 14.4. Detection limit

#### Semi-quantitative protocol

**Calibrator measurement range:** 12.5% - 100 %

**Limit of Detection (LoD) = 8%**

## 15. TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSES	SOLUTION
<b>Control values out of range</b>	Incorrect temperature, timing or pipetting, reagents are not mixed	Check that the time and temperature were correct. Repeat test.
	Cross contamination of controls	Pipette carefully.
	Optical pathway is not clean.	Check for dirt or air-bubbles in the wells. Wipe plate bottom and reread.
	Controls (positive and/or activity controls) are not correctly reconstituted. Improper dilution of calibrator.	Check the controls, dissolve a new. Check the preparation and make a new dilution.
<b>All test results negative</b>	One or more reagents are not added or added in wrong sequence.	Recheck procedure. Check for unused reagents. Repeat test.
	Antigen coated plate is inactive.	Check for obvious moisture in unused wells. Wipe plate bottom and reread.
	Serum inactive.	Dilute new samples.
<b>All sample wells are visibly yellow.</b>	Contaminated buffers or reagents.	Check all solutions for turbidity.
	Washing solution is contaminated.	Use clean container. Check the quality of water used for preparation of solution.
	Improper dilution of serum.	Repeat test.
<b>Poor precision.</b>	Pipette delivery CV >5% or samples not mixed.	Check the calibration of pipette. Use reproducible technique. Avoid air bubbles in pipette tip.
	Serum or reagents are not mixed sufficiently or not equilibrated to room temperature.	Mix all reagents gently but thoroughly and equilibrate to room temperature.
	Reagent addition is taking too long time, inconsistency in timing intervals.	Develop consistent uniform technique and use multi-tip device or auto-dispenser to decrease time.
	Optical pathway not clean.	Check for air-bubbles in the wells. Wipe plate bottom and reread.
	Washing not consistent, trapped bubbles, washing solution left in the wells.	Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in the well. After last wash, empty the wells by tapping the strip on an absorbent tissue.

## 16. LITERATURE REFERENCES

- Walport M. Complement (First of two parts). *N Engl J Med* 2001; 344:1058-66.
- Walport M. Complement (Second of two parts). *N Engl J Med* 2001; 344:1140-44.
- Roos A *et al.* Functional characterization of the lectin pathway of complement in human serum. *Mol Immunol* 2003; 39:655-68.
- Fredrikson GN *et al.* New procedure for the detection of complement deficiency by ELISA. Analysis of activation pathways and circumvention of rheumatoid factor influence. *J Immunol Methods*. 1993 Dec 3;166(2):263-70.
- Seelen MA *et al.* Functional analysis of the classical, alternative and MBL pathways of the complement system: standardization and validation of a simple ELISA. *J Immunol Meth* 2005; 296:187-98

- Salvador-Morales C, Sim RB. Handbook of Immunological Properties of Engineered Nanomaterials. 2013, 1st Ed, World Scientific Publishing (ISBN: 978-4390-25-5).
- Tudoran R & Kirschfink M. Modern Complement analysis: indications, methods and outlook. *LaboratoriumsMedizin* 2012; 36(3):--.
- Botto M *et al.* Complement in human disease: Lessons from complement deficiencies. *Mol Immunol* 2009; 46:2774-83.
- Mollnes TE *et al.* Complement analysis in the 21st century. *Mol Immunol* 2007; 44:3838-49.
- Nilsson B, Nilsson Ekdahl K. Complement Diagnostics: Concepts, Indications, and Practical Guidelines. *Clin Develop Immunol*; 2012, Art ID 962702.

## 17. DESCRIPTION OF SYMBOLS

The following symbols may appear on the packaging and the label:

	Batch code
	Catalog number.
	Use-by-date.
	Temperature limit.
	Biological risk.
	Consult instructions for use.
	In vitro diagnostic use.
	Manufacturer.
	Contains sufficient for 96 tests.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive.
	Health hazard

Ag	Antigen (coated plate).
DIL	Diluent.
BUF WASH 30X	Wash buffer, 30x concentrate.
CONTROL + LYO	Positive control.
CONTROL AC LYO	Lyophilized activity control.
CONTROL -	Negative control.
CONJ	Conjugate.
SUBS pNPP	Substrate pNPP.
CH REP	Swiss Representative.

## 18. DOCUMENT HISTORY

Version	Changes
LABEL-DOC-0031 v 6.0	<p>Major editorial changes with new document layout. The procedures for semi-quantitative and qualitative application have been merged and separate sub-headings have been introduced where applicable.</p> <p>Information that “The analysis should be performed by trained laboratory professionals” added.</p> <p>Information regarding metrological traceability added.</p> <p>Specification for positive control is changed from &gt;1 to &gt;1.0 for clarification.</p> <p>Information regarding CO<sub>2</sub> cabinet added.</p>
LABEL-DOC-0031 v 7.0	<p>Increased transparency in section 14.1 that the semi-quantitative data is based on mathematical theoretical calculation.</p> <p>Values in table 6 and 10 were updated to correspond to the performance evaluation published in version 5.0.</p>

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## 16. ODKAZY NA LITERATURU

- Walport M. Complement (First of two parts). N Engl J Med 2001; 344:1058-66.
- Walport M. Complement (Second of two parts). N Engl J Med 2001; 344:1140-44.
- Roos A *et al.* Functional characterization of the lectin pathway of complement in human serum. Mol Immunol 2003; 39:655-68.
- Fredrikson GN *et al.* New procedure for the detection of complement deficiency by ELISA. Analysis of activation pathways and circumvention of rheumatoid factor influence. J Immunol Methods. 1993 Dec 3;166(2):263-70.
- Seelen MA *et al.* Functional analysis of the classical, alternative and MBL pathways of the complement system: standardization and validation of a simple ELISA. J Immunol Meth 2005; 296:187-98

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- Walport M. Complement (First of two parts). N Engl J Med 2001; 344:1058-66.
- Walport M. Complement (Second of two parts). N Engl J Med 2001; 344:1140-44.
- Roos A *et al.* Functional characterization of the lectin pathway of complement in human serum. Mol Immunol 2003; 39:655-68.
- Fredrikson GN *et al.* New procedure for the detection of complement deficiency by ELISA. Analysis of activation pathways and circumvention of rheumatoid factor influence. J Immunol Methods. 1993 Dec 3;166(2):263-70.
- Seelen MA *et al.* Functional analysis of the classical, alternative and MBL pathways of the complement system: standardization and validation of a simple ELISA. J Immunol Meth 2005; 296:187-98
- Salvador-Morales C, Sim RB. Handbook of Immunological Properties of Engineered Nanomaterials. 2013, 1st Ed, World Scientific Publishing (ISBN: 978-4390-25-5).
- Tudoran R & Kirschfink M. Modern Complement analysis: indications, methods and outlook. LaboratoriumsMedizin 2012; 36(3):--.
- Botto M *et al.* Complement in human disease: Lessons from complement deficiencies. Mol Immunol 2009; 46:2774-83.
- Mollnes TE *et al.* Complement analysis in the 21st century. Mol Immunol 2007; 44:3838-49.
- Nilsson B, Nilsson Ekdahl K. Complement Diagnostics: Concepts, Indications, and Practical Guidelines. Clin Develop Immunol; 2012, Art ID 962702.

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- Walport M. Complement (First of two parts). N Engl J Med 2001; 344:1058-66.
- Walport M. Complement (Second of two parts). N Engl J Med 2001; 344:1140-44.
- Roos A *et al.* Functional characterization of the lectin pathway of complement in human serum. Mol Immunol 2003; 39:655-68.
- Fredrikson GN *et al.* New procedure for the detection of complement deficiency by ELISA. Analysis of activation pathways and circumvention of rheumatoid factor influence. J Immunol Methods. 1993 Dec 3;166(2):263-70.
- Seelen MA *et al.* Functional analysis of the classical, alternative and MBL pathways of the complement system: standardization and validation of a simple ELISA. J Immunol Meth 2005; 296:187-98
  
- Salvador-Morales C, Sim RB. Handbook of Immunological Properties of Engineered Nanomaterials. 2013, 1st Ed, World Scientific Publishing (ISBN: 978-4390-25-5).
- Tudoran R & Kirschfink M. Modern Complement analysis: indications, methods and outlook. LaboratoriumsMedizin 2012; 36(3):--.
- Botto M *et al.* Complement in human disease: Lessons from complement deficiencies. Mol Immunol 2009; 46:2774-83.
- Mollnes TE *et al.* Complement analysis in the 21st century. Mol Immunol 2007; 44:3838-49.
- Nilsson B, Nilsson Ekdahl K. Complement Diagnostics: Concepts, Indications, and Practical Guidelines. Clin Develop Immunol; 2012, Art ID 962702.

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- Walport M. Complement (First of two parts). *N Engl J Med* 2001; 344:1058-66.
- Walport M. Complement (Second of two parts). *N Engl J Med* 2001; 344:1140-44.
- Roos A *et al.* Functional characterization of the lectin pathway of complement in human serum. *Mol Immunol* 2003; 39:655-68.
- Fredrikson GN *et al.* New procedure for the detection of complement deficiency by ELISA. Analysis of activation pathways and circumvention of rheumatoid factor influence. *J Immunol Methods*. 1993 Dec 3;166(2):263-70.
- Seelen MA *et al.* Functional analysis of the classical, alternative and MBL pathways of the complement system: standardization and validation of a simple ELISA. *J Immunol Meth* 2005; 296:187-98
- Salvador-Morales C, Sim RB. *Handbook of Immunological Properties of Engineered Nanomaterials*. 2013, 1st Ed, World Scientific Publishing (ISBN: 978-4390-25-5).
- Tudoran R & Kirschfink M. Modern Complement analysis: indications, methods and outlook. *LaboratoriumsMedizin* 2012; 36(3):--.
- Botto M *et al.* Complement in human disease: Lessons from complement deficiencies. *Mol Immunol* 2009; 46:2774-83.
- Mollnes TE *et al.* Complement analysis in the 21st century. *Mol Immunol* 2007; 44:3838-49.
- Nilsson B, Nilsson Ekdahl K. Complement Diagnostics: Concepts, Indications, and Practical Guidelines. *Clin Develop Immunol*; 2012, Art ID 962702.

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
### 3. PRINCIP WIESLAB® COMPLEMENT SYSTEM CLASSICAL PATHWAY

Test WIESLAB® Complement system Classical pathway kombinuje principe hemolitičkog testa za aktivaciju komplementa uz upotrebu označenih antitela specifičnih za neoantigen proizveden kao rezultat aktivacije komplementa. Količina nastalog neoantigena je proporcionalna funkcionalnoj aktivnosti klasičnog puta.

U kompletu klasičnog puta sistema komplementa, bunarčići traka za mikrotitraciju su obloženi specifičnim aktivatorom klasičnog puta. Ploča u kombinaciji sa sastavom pufera za razblaživanje uzorka i nivoom razblaživanja seruma pacijenta osigurava aktiviranje samo klasičnog puta. Tokom inkubacije razblaženog seruma pacijenta u bunarčićima, komplement se aktivira specifičnim premazom.

Bunarčići se zatim ispiraju i količina kompleksa C5b-9 formirana na površini ploče se detektuje antitelom označenim alkalnom fosfatazom specifičnim za neoantigen C5b-9 formiran tokom formiranja membranolitičkog kompleksa (MAC).

Nakon daljeg koraka pranja, detekcija specifičnih antitela se dobija inkubacijom sa rastvorom supstrata alkalne fosfataze. Količina aktivacije komplementa u vezi je sa intenzitetom boje i meri se u smislu apsorbance (optička gustina, OD).

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**Inspira GmbH**
**SVAR LIFE SCIENCE AB**