

Certificate of Analysis

IMMUNOSCAN CCPlus®

REF RA-96 PLUS

LOT XS 1132


2019-08-12

Reagents			LOT		Specification	Result
Ag			US 2936	2019-11-13		
CAL	A		XS 1133	2019-08-12	OD	1.76
CAL	B		XS 1134	2019-08-12	OD	1.48
CAL	C		XS 1135	2019-08-12	OD	0.90
CAL	D		XS 1136	2019-08-12	OD	0.36
CAL	E		XS 1137	2019-08-12	OD ≥0.15	0.19
CONTROL	REF		XS 1138	2019-08-12	OD ≥0.15	0.19
CONTROL	+		XS 1139	2019-08-12	180-340 U/mL	252
					Ratio 4.2-7.6	5.3
CONTROL	-		XS 1140	2020-01-28	OD	0.01
					Ratio <0.95	0.05
CONJ			XS 1141	2019-08-12		
DIL			XS 1143	2020-01-28	OD	0.05
BUF	WASH	20X	XS 1142	2020-01-29		
SUBS		TMB	XS 1058	2020-01-31		
H ₂ SO ₄		0.5M	XS 1200	2020-01-23		

It is hereby certified that this LOT has been manufactured, tested and approved according to our procedures and in compliance with ISO 13485.

Date of approval: 2018-04-24



Quality Control



Quality Assurance

