

Certificate of Analysis

IMMUNOSCAN CCPlus®

REF RA-96 PLUS

LOT XS 1730


2019-11-01

Reagents			LOT		Specification	Result
Ag			XS 1393	2020-03-10		
CAL	A		XS 1731	2019-11-04	OD	1.82
CAL	B		XS 1732	2019-11-04	OD	1.51
CAL	C		XS 1733	2019-11-04	OD	0.92
CAL	D		XS 1734	2019-11-04	OD	0.36
CAL	E		XS 1735	2019-11-04	OD ≥ 0.15	0.20
CONTROL	REF		XS 1736	2019-11-04	OD ≥ 0.15	0.19
CONTROL	+		XS 1737	2019-11-04	180-340 U/mL	249
					Ratio 4.2-7.6	5.3
CONTROL	-		XS 1738	2020-05-01	OD	0.01
					Ratio < 0.95	0.05
CONJ			XS 1740	2019-11-01		
DIL			XS 1729	2020-05-01	OD	0.05
BUF	WASH	20X	XS 1728	2020-04-29		
SUBS		TMB	XS 1741	2020-04-30		
H ₂ SO ₄		0.5M	XS 1739	2020-04-29		

It is hereby certified that this LOT has been manufactured, tested and approved according to our procedures and in compliance with ISO 13485.

Date of approval: 2018-06-14

Gordana Prijic Priban
Quality Control

Belma Tinjak
Quality Assurance

