

Certificate of Analysis

IMMUNOSCAN CCPlus®

REF RA-96 PLUS

LOT XS 2590  2019-09-16

Reagents		LOT			Specification	Result
Ag		XS 1393	2020-03-10			
CAL	A	XS 1731	2019-11-04		OD	1.82
CAL	B	XS 1732	2019-11-04		OD	1.51
CAL	C	XS 1733	2019-11-04		OD	0.92
CAL	D	XS 1734	2019-11-04		OD	0.36
CAL	E	XS 1735	2019-11-04		OD ≥0.15	0.20
CONTROL	REF	XS 1736	2019-11-04		OD ≥0.15	0.19
CONTROL	+	XS 1737	2019-11-04		180-340 U/mL	249
CONTROL	-	XS 1738	2020-05-01		Ratio 4.2-7.6	5.3
					OD	0.01
					Ratio <0.95	0.05
CONJ		XS 1740	2019-11-01			
DIL		XS 1729	2020-05-01		OD	0.05
BUF	WASH	20X	XS 1728	2020-04-29		
SUBS	TMB	XS 1741	2020-04-30			
H ₂ SO ₄	0.5M	XS 1739	2020-04-29			

It is hereby certified that this LOT has been manufactured, tested and approved according to our procedures and in compliance with ISO 13485.

Date of approval: 2018-09-26

Anna Nilsen

Quality Control

Netta Olofsson

Quality Assurance

